

# Appendix A1

January 23, 2023

## 1 Appendix A.1 - Data Collection

```
[ ]: from pydaqmx_helper.adc import ADC
     from pydaqmx_helper.digital_io import Digital_IO
     import time
```

```
[ ]: import matplotlib.pyplot as plt
     import numpy as np
```

```
[ ]: adc = ADC()
```

```
[ ]: adc.addChannels([1])
     print(adc.getActiveChannels())
```

Activated Channel 1  
[1]

```
[ ]: dio = Digital_IO('0', 'output')
```

Created digital output port: Dev2/port0

```
[ ]: def SetState(state, dio=dio):
     dio.write(state)

     def Turn(enabled, N=1, prism=False, detector=False, direction="cw",
             ↪doubleDetector=False, delay=0.1, \
             reverse=False, stationary=False):

         direction = direction.lower()
         voltages = []*N
         uncertainties = []*N
         iterations = []*N

         # Enable
         if not stationary:
             if enabled:
                 currentState = 1
```

```

        SetState(currentState)
    else:
        currentState = 1
        SetState(0)
        return
    elif stationary:
        currentState = 0
        SetState(currentState)

    # Selecting Direction
    if (direction == "ccw"):
        currentState += 2
        SetState(currentState)

    elif (direction == "cw"):
        SetState(currentState)

    else:
        currentState = 0
        SetState(currentState)
        print("Motor Disabled, invalid input")
        return

    print("Starting Cycle...")
    # Do turn cycle
    # Delay 100 ms
    time.sleep(0.1)

    for i in range(N):
        print(f"Cycle {i}/{N}", end="\r")

        if detector:
            # Turn motor two
            currentState += 16
            SetState(currentState)
            currentState -= 16
            SetState(currentState)

            time.sleep(delay)

        if doubleDetector:
            currentState += 16
            SetState(currentState)
            currentState -= 16
            SetState(currentState)
            time.sleep(delay)

```

```

        if prism:
            # Turn motor one
            currentState += 8
            SetState(currentState)
            currentState -= 8
            SetState(currentState)

    if prism and not detector:
        # Turn motor one
        currentState += 8
        SetState(currentState)
        currentState -= 8
        SetState(currentState)

    time.sleep(delay)

time.sleep(delay)

iterations.append(i)
V=adc.sampleVoltages(100, 10000)[1] # 100 samples at 10000 Hz
voltages.append(np.mean(V))
uncertainties.append(np.std(V))

if stationary:
    time.sleep(1)

print("Cycle Complete, disabling motors")
currentState = 0
SetState(currentState)

if reverse:
    reverseIterations = iterations[::-1]
    return [reverseIterations, voltages, uncertainties]

return [iterations, voltages, uncertainties]

```

```
[ ]: Turn(False)
```

```
[ ]: Turn(True, 10, prism=True, detector=False, direction="cw", doubleDetector=True);
```

Starting Cycle...

Cycle Complete, disabling motors

```
[ ]: # Calibrate turns to degrees function
      # From calibration, the system moves 45 degrees in 1000 turns for both the
      ↪prism and the detector
```

```
def DegToTurns(deg):
    ratio = 45/1000 # degrees / turns
    return deg / ratio

def TurnsToDeg(turns):
    ratio = 45/1000 # degrees / turns
    return turns * ratio
```

```
[ ]: DegToTurns(65)
```

```
[ ]: 1444.4444444444446
```

```
[ ]: print(f"Red and Green Offset: {TurnsToDeg(670)}")
print(f"Blue Offset: {TurnsToDeg(970)}")
```

Red and Green Offset: 30.15

Blue Offset: 43.65

```
[ ]: def WriteFile(fileName, data):

    data = np.transpose(data)

    with open(fileName, "w") as f:
        np.savetxt(fileName, data)
```

```
[ ]: def LoadFile(path):
    data = np.array(np.loadtxt(path))
    return data
```

```
[ ]: def PlotFile(path):
    data = LoadFile(path)

    plt.xlabel("N")
    plt.ylabel("Intensity")
    plt.grid()
    plt.plot(data[:,0], data[:,1])
```

## 2 Data Acquisition

```
[ ]: # Note that we assume detector needs to move twice as much as the prism

# 0 deg is defined at perfect reflection, this is read at an angle of 130 deg
# Starting with the prism at 120
```

### 2.0.1 Example data acquisition

```
[ ]: WriteFile("15nm/red/red01.txt", Turn(True, 600, prism=True, detector=True, ↵  
↵direction="cw", doubleDetector=True))
```

Starting Cycle...

Cycle Complete, disabling motors

```
[ ]: PlotFile("15nm/red/red01.txt")
```

