### **CCDs**

#### October 23, 2022

```
[1]: # IMPORTS #

# For handeling FITS files
from astropy.io import fits
import scipy
from glob import glob

import numpy as np
from matplotlib import pyplot as plt

[2]: # The below commands make the font and image size bigger
plt.rcParams.update({'font.size': 26})
plt.rcParams["figure.figsize"] = (12,12)
```

# 1 Determining Readout Noise and Gain

Τ

```
BITPIX =
                            16
NAXIS
                             2
NAXIS1 =
                           100
NAXIS2 =
                           100
BZERO =
                    32768.000000
BSCALE =
                      1.000000
                      0.000000
DATAMIN =
DATAMAX =
                    65535.000000
INSTRUME=
                    'ATIK-314L: fw rev 3.01'
FILTER =
                    'Position 1'
EXPTIME =
                    0.001
                    '2022-10-05T14:01:39.2'
DATE-OBS=
XPIXSZ =
                    6.450
                    6.450
YPIXSZ =
XBINNING=
YBINNING=
```

SIMPLE =

```
XORGSUBF=
                         661
    YORGSUBF=
                         489
    XPOSSUBF=
                        661
    YPOSSUBF=
                        489
    CBLACK =
                        216
    CWHITE =
                        297
    CCD-TEMP=
                        -7.5
    SWCREATE=
                         'Artemis Capture'
    END
[4]: fIm = fits.open(r"/home/daraghhollman/jupyter/UCD_PASS_Labs/CCDs/Data_D1/

daragh_ccd/flat/flat_03.fit")
     print(fIm[0].header)
    SIMPLE =
                                  Τ
    BITPIX =
                                 16
                                  2
    NAXIS
    NAXIS1 =
                                100
    NAXIS2 =
                                100
    BZERO =
                        32768.000000
    BSCALE =
                           1.000000
    DATAMIN =
                           0.000000
    DATAMAX =
                        65535.000000
    INSTRUME=
                         'ATIK-314L: fw rev 3.01'
    FILTER =
                         'Position 1'
    EXPTIME =
                        0.100
    DATE-OBS=
                         '2022-10-05T14:29:49.7'
    XPIXSZ =
                        6.450
    YPIXSZ =
                        6.450
                         1
    XBINNING=
    YBINNING=
                         1
    XORGSUBF=
                        661
    YORGSUBF=
                        489
    XPOSSUBF=
                        661
    YPOSSUBF=
                        489
    CBLACK =
                        16450
    CWHITE =
                         18030
    CCD-TEMP=
                        -7.2
    SWCREATE=
                         'Artemis Capture'
    END
[5]: biasFilesA = glob(r"/home/daraghhollman/jupyter/UCD_PASS_Labs/CCDs/Data_D1/

¬daragh_ccd/bias/bias*")
     biasFilesB = glob(r"/home/daraghhollman/jupyter/UCD_PASS_Labs/CCDs/Data_D1/

daragh_ccd/bias_1/bias*")

     biasFITSA = [fits.open(el) for el in biasFilesA] # A list of 30 fits files
```

```
biasFITSB = [fits.open(el) for el in biasFilesB] # A list of 30 fits files
biasFITS = biasFITSA+biasFITSB
```

```
[8]: # Biases
     i=0
     myFlats = list()
     while i < len(flatFITS):</pre>
         # Check if there are still pairs left
         if (len(flatFITS) - i) >= 2:
             myFlats.append([flatFITS[i], flatFITS[i+1]])
             i+=2
     j = 0
     myBiases = list()
     while j < len(biasFITS):</pre>
         # Check there still exist pairs
         if (len(biasFITS) - i >= 2):
             myBiases.append([biasFITS[j], biasFITS[j+1]])
             j+=2
     def gainAndReadnoise(flats, biases, std=False):
         # where flats and biases are lists with elements containing pairs of flats,
      ⇔or biases
         gainList = list()
```


```
readnoiseList = list()
   for flatPair, biasPair in zip(flats, biases):
        meanFlats = [flatPair[0][0].data.mean(), flatPair[1][0].data.mean()]
        meanBiases = [biasPair[0][0].data.mean(), biasPair[1][0].data.mean()]
        flatDiff = abs(flatPair[0][0].data - flatPair[1][0].data)
        biasDiff = abs(biasPair[0][0].data - biasPair[1][0].data)
        std_flatDiff = flatDiff.std()
        std_biasDiff = biasDiff.std()
       gain = ((meanFlats[0] + meanFlats[1]) - (meanBiases[0] +__
 meanBiases[1])) / (std_flatDiff**2 - std_biasDiff**2)
        readnoise = gain * std_biasDiff / np.sqrt(2)
        gainList.append(gain)
        readnoiseList.append(readnoise)
   averageGain = np.mean(gainList)
    std gain = np.std(gainList)
   averageReadnoise = np.mean(readnoiseList)
    std_readnoise = np.std(readnoiseList)
   if std==True: return [[averageGain, averageReadnoise], [std_gain,_
 →std_readnoise]]
    else: return [averageGain, averageReadnoise]
res = gainAndReadnoise(myFlats, myBiases, std=True)
gain = res[0][0]
readnoise = res[0][1]
print(f"Gain: {gain} with standard deviation: {res[1][0]}")
print(f"Readnoise: {readnoise} with standarad deviation: {res[1][1]}")
```

Gain: -0.001577846989986806 with standarad deviation: 0.000748331061405 Readnoise: -36.53031784248212 with standarad deviation: 17.324857482162162

# 2 Spectrum of the Arc Lamp

```
plt.ylabel("y coordinate")
```

## [9]: Text(0, 0.5, 'y coordinate')




```
[10]: line_values = arcImage[0].data[800]
print(max(line_values))
```

#### 52351


```
[12]: def gaussian(x, mu, std, scale):
          ans = list()
          for el in x:
              ans.append(scale * np.exp(-(el-mu)**2/(2*std**2)))
          return ans
      pars01, cov01 = scipy.optimize.curve_fit(gaussian, x[0:150], y[0:150], [80, 10, ___
       →10000])
      pars02, cov02 = scipy.optimize.curve_fit(gaussian, x[260:290], y[260:290],__
       \hookrightarrow [278, 4, 52000])
      pars03, cov03 = scipy.optimize.curve fit(gaussian, x[920:980], y[920:980],
       \hookrightarrow [950, 4, 12000])
      pars04, cov04 = scipy.optimize.curve_fit(gaussian, x[1050:1200], y[1050:1200],
       \hookrightarrow[1125, 10, 400])
      meanX = [pars01[0], pars02[0], pars03[0], pars04[0]]
[13]: # First peak : [0:150]
      # Second peak : [200:350]
      # Third peak : [800:1100]
      # Fourth peak : [1050:1200]
      plt.plot(x, y, color="black", label="Data")
      # Plan is to plot vertical lines displaying the estimated x and the uncertainty
      # Peak 1
      plt.plot([65]*11000, np.arange(0, 11000, 1), color="red")
      # Peak 2
      plt.plot([277]*56000, np.arange(0, 56000, 1), color="red", label="Line peak")
      plt.plot([965]*15000, np.arange(0, 15000, 1), color="red")
      plt.plot([1145]*6000, np.arange(0, 6000, 1), color="red")
      plt.legend()
      plt.grid()
      plt.xlabel("x coordinate of pixel")
      plt.ylabel("Counts")
```

[13]: Text(0, 0.5, 'Counts')



```
plt.legend()
plt.grid()
plt.ylabel("x pixel coordinate")
plt.xlabel("wavelength (nm)")
```

[14]: Text(0.5, 0, 'wavelength (nm)')




# 3 Investigate the Resolution of the Spectrograph

```
(1039, 1391)
```

```
[16]: plt.imshow(imageNa[0].data, cmap="Greys_r")
plt.plot(np.arange(0, 1391, 1), [800]*(1391), color="red")
plt.colorbar()
```

[16]: <matplotlib.colorbar.Colorbar at 0x7f90a4fb7be0>



```
[17]: line_valuesNa = imageNa[0].data[800]
print(max(line_valuesNa))
```

45180

```
[18]: xNa = list()
yNa = list()
k = 0
```


```
for el in line_valuesNa:
    xNa.append(k)
    yNa.append(el)
    k+=1
```

```
[19]: plt.plot(xNa, yNa, color="black", label="Data")

# Peak 1
plt.plot([120]*48000, np.arange(0, 48000, 1), color="red")
# Peak 2
plt.plot([465]*8000, np.arange(0, 8000, 1), color="red", label="Line peak")
# Peak 3
plt.plot([575]*10000, np.arange(0, 10000, 1), color="red")
# Peak 4
plt.plot([1045]*6000, np.arange(0, 6000, 1), color="red")
# Peak 5
plt.plot([1070]*5000, np.arange(0, 5000, 1), color="red")


plt.legend()
plt.grid()
plt.xlabel("x coordinate of pixel")
plt.ylabel("Counts")
```

[19]: Text(0, 0.5, 'Counts')



```
A = C / (abs(std) * np.sqrt(2 * np.pi))
   B = np.exp(-0.5 * ((x - mean) / abs(std))**2)
   return A*B + H
gaussPars01, gaussCov01 = scipy.optimize.curve_fit(gaussian, NaWavelengths[100:
 →150], yNa[100:150], \
                                          [565, 10, 40000, 5000])
gaussPars02, gaussCov02 = scipy.optimize.curve_fit(gaussian, NaWavelengths[550:
 →600], yNa[550:600], \
                                          [495, 5, 8000, 5000])
gaussPars03, gaussCov03 = scipy.optimize.curve fit(gaussian, NaWavelengths[1020:
 →1065], yNa[1020:1065], \
                                          [410, 5, 4000, 3000])
gaussPars04, gaussCov04 = scipy.optimize.curve fit(gaussian, NaWavelengths[450:
 →480], yNa[450:480], \
                                          [520, 5, 6000, 5000])
plt.plot(NaWavelengths[100:150], gaussian(NaWavelengths[100:150],
 ⇒gaussPars01[0], gaussPars01[1], gaussPars01[2], \
                                     gaussPars01[3]), label=f"Gaussian atu
plt.plot(NaWavelengths[550:600], gaussian(NaWavelengths[550:600],
 ⇒gaussPars02[0], gaussPars02[1], gaussPars02[2],
                                     gaussPars02[3]), label=f"Gaussian at_
 plt.plot(NaWavelengths[1020:1065], gaussian(NaWavelengths[1020:1065],
 ⇒gaussPars03[0], gaussPars03[1], gaussPars03[2], \
                                     gaussPars03[3]), label=f"Gaussian at⊔
 plt.plot(NaWavelengths[450:480], gaussian(NaWavelengths[450:480],
 →gaussPars04[0], gaussPars04[1], gaussPars04[2], \
                                     gaussPars04[3]), label=f"Gaussian atu
 plt.legend()
plt.grid()
plt.xlabel("Wavelength (nm)")
plt.ylabel("Counts")
```

[20]: Text(0, 0.5, 'Counts')



#### 3.1 Find Resolution

```
[21]: def resolution(std, wavelength):
    dL = 2.355 * std
    res = wavelength / dL
    return res


[22]: def resUncertainty(wavelength, FWHM, uW, uFWHM):
    return np.sqrt((1/FWHM)**2 * uW**2 + (wavelength/(FWHM**2))**2 * uFWHM**2)

[23]: point1 = resolution(gaussPars01[1], gaussPars01[0])
    point2 = resolution(abs(gaussPars04[1]), gaussPars04[0])
    point3 = resolution(gaussPars02[1], gaussPars02[0])
    point4 = resolution(gaussPars03[1], gaussPars03[0])
    resolutions = [point1, point2, point3, point4]
```

Resolution of the spectrograph at 571.39 nm: 301.36518 Resolution of the spectrograph at 514.48 nm: 284.79700 Resolution of the spectrograph at 496.73 nm: 304.49310 Resolution of the spectrograph at 419.06 nm: 206.63528


[17.86805221558772, 34.89525500297671, 38.51343788860465, 44.725515255856266]

- 0.628751594086656 -40.313700669343525
- 0.2629224814699645 132.34737527398428



## 4 Effect of slit width on the resolution


[27]: <matplotlib.colorbar.Colorbar at 0x7f90a493fe80>



```
[28]: line_valuesSlit = imageSlit[0].data[800]
print(max(line_valuesSlit))
```

#### 24947

### [30]: Text(0, 0.5, 'Counts')



Note that the peaks are all wider which will yield a larger standard deviation and hence a smaller resolution, accuracy decreases. The images were taken at different exposures to avoid peaking the data however it is clear that the larger slit has wider peaks.

[]:[